

New Product Release

'Ready to Use' Ceramic Crucibles



We are pleased to announce a new ceramic crucible for inorganic Carbon and Sulphur analysis for metals and other inorganic materials. The crucible is compatible with induction furnace instruments from Leco[®], Eltra[®], Horiba[®], Bruker[®] and others.

This crucible is designed to be used straight from the box without requiring pre-burning or other pre-treatment in most applications saving substantial time and cost. The product is supplied foil-wrapped in boxes of 500. A sample pack of 20 crucibles is also available, ask our sales team for more information.

Typical crucible blank values are less than 10ppm Carbon, and less than 5ppm Sulphur.

| Part Number | Description | Pack quantity |
|-------------|--|---------------|
| C4516 | Ready to use ceramic crucible | 500 |
| Z2038 | Sample Pack of C4516 Ready to use ceramic crucible | 20 |



These products offer the same high quality and value for money that is typical of Elemental Microanalysis' entire product range.

Available to order from our ISO 9001 factory [online](#), [by email](#), telephone or fax. Elemental Microanalysis offers extensive knowledge of analysis, over 40 years of experience and a worldwide reputation for excellence.