

## New Product Range

### Product range suitable for ELTRA® Elementrac ON/OH/OHN-p


Elemental Microanalysis is pleased to introduce a new range of products suitable for the ELTRA® Elementrac ON/OH/OHN-p analyser.

Renowned for longevity of service, the robust Elementrac ON/OH/OHN -p analyser from ELTRA® measures oxygen, nitrogen and hydrogen concentrations in inorganic samples like steel, iron, copper or ceramics.

The high quality compatible products we supply to support the continued running of your Elementrac analyser are produced in our ISO 9001 facility, we guarantee that our products are as good or better than those from your current supplier.

#### Reagent and catalyst tubes

ELTRA® part	EMAL part number	Reagent and catalyst tubes
20040	C6506	Catalyst tube 200mm x 17mm
11064 88400-0006	C6501	Reagent tube 280mm x 16mm
11064-3001	C6521	Reagent tube 120mm x 16mm ID
88400-0422	C6508	Reagent tube



#### Operational tools

ELTRA® part	EMAL part number	Operational tools
31360	E2421	Graphite tip
31365	E2420	Graphite tip holder
71029	E2423	Cleaning brush for furnace
71032	E2755	Brush
92610	E2158	Silicone grease 100gm
88400-0229	E2472	Crucible tweezers
	E2480	Cleaning brush assembly
	E2756	Brush handle



Graphite tip  
E2421



Graphite tip holder E2420



Crucible  
tweezers  
E2472



Cleaning brush E2423



Brush E2755

The use of ELTRA® part numbers is intended for convenience only, and does not imply the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All ELTRA® trademarks acknowledged.

Product range suitable for ELTRA®  
Elementrac ON/OH/OHN-p cont'd

## Graphite crucibles



Outer graphite crucible C4610, C4615



Inner graphite crucible C4611, C4609, C4614



Graphite crucible C4612, C4637, C4638



Inner graphite crucible C4611, C4609, C4614

ELTRA® part	EMAL part number	Graphite crucibles	Pack
90180	C4611	Inner graphite crucible for use with C4610	100
90185	C4610	Outer graphite crucible for use with C4609, C4611	50
90185	C4615	Outer graphite crucible for use with C4609, C4611	500
90190	C4612	Graphite crucible	100
	C4609	Inner graphite crucible for use with C4610	100
	C4614	Inner graphite crucible for use with C4610	1000
	C4637	Graphite crucible	1000
	C4638	Graphite crucible	5000

## O-rings

ELTRA® part	EMAL part number	O-rings
70210	E1214	O-ring
70230	E1151	O-ring 9x3
70410	E1157	O-ring



E1214



E1151



E1157

The use of ELTRA® part numbers is intended for convenience only, and does not imply the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All ELTRA® trademarks acknowledged.



Product range suitable for ELTRA®  
Elementrac ON/OH/OHN-p cont'd

## Reagents and fillings

ELTRA® part	EMAL part number	Reagents and fillings	Pack
90200	B1161	EMADrone magnesium perchlorate	454g
90210	B1164	EMASorb B granular 1.6 to 3mm 7 to 14 mesh	500g
90210	B1317	EMASORB B II (8-20mesh)	500g
90270	B1519	Schutze reagent	100g
90289	B1123	Rare earth copper oxide granular 0.85 to 1.7mm	50g
90330	B1350	Silica (quartz) wool fine	50g
90331	B1166	Glass wool borosilicate fibres	454g
90331	B1266	Glass wool rope borosilicate fibres	454g
90800	B1516	Graphite powder	100g
	B1162	EMASorb A granular 0.8 to 1.6mm 14 to 25 mesh	500g
	B1167	Glass wool borosilicate fibres	227g
	B1267	Glass wool rope borosilicate fibres	227g
	B1316	EMASORB A II (20-30 mesh)	500g



EMADrone magnesium perchlorate B1161



EMASorb B granular B1164



EMASORB B II B1317



Rare earth copper oxide B1123



Schutze reagent B1519

These products offer the same high quality and value for money that is typical of Elemental Microanalysis' entire product range.

Available to order from our [ISO 9001](#) factory online, by email, telephone or fax. Elemental Microanalysis offers extensive knowledge of analysis, over 40 years of experience and a worldwide reputation for excellence.

*The use of ELTRA® part numbers is intended for convenience only, and does not imply the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All ELTRA® trademarks acknowledged.*

### Elemental Microanalysis Ltd

1 Hameldown Road, Okehampton, EX20 1UB, UK

t. +44(0)1837 54446 • f. +44(0)1837 54544

[www.elementalmicroanalysis.com](http://www.elementalmicroanalysis.com) • [enquiries@microanalysis.co.uk](mailto:enquiries@microanalysis.co.uk)



Available from your local Dealer:

