

New Product Information

New products for EuroVector EA3100 Plus

Elemental Microanalysis and EuroVector are pleased to introduce a range of new products for the EA3100 Plus analyser, featuring a high thickness reactor tube with 3 mm walls, C1151V, that offers extra durability and higher throughput of samples.

The new reactor tube has a larger diameter (25 mm) than the reactor tube used in the EA3100, but users of the EA3100 can upgrade to the 25 mm configuration, with some adjustment first by a suitably qualified engineer.

The comprehensive new product range also includes prepacked reactors for CHNS analysis (X1076V), or oxygen pyrolysis (X1079V), two start-up consumables kits, and two consumables kits for continuing analysis. The full contents of these kits is shown below.



C1151V
high thickness
reactor tube

EMAL part	EuroVector part	Description
C1151V	C09-013	Quartz empty reactor tube high thickness 25/19 mm
X1076V	C10-151	Packed reactor 25/19 mm CHNS (1pc)
X1079V	C10-205	Prepacked oxygen pyrolysis tube 25/19 mm (1pc)
X3276	C11-078	Start up kit for 1000 analyses 25/19 mm (CHNS w heteroatoms): <ul style="list-style-type: none"> X1076V (1 each): Packed reactor 25/19 mm CHNS (1pcs) C10-151 D1008 (4 each): Tin capsules pressed standard weight 8 x 5 mm, pack of 250 B2131 (1 each): Standards kit CHNS/CHNS-O, cystine, sulphanilamide, BBOT, sulphanilic acid, 4 x 0.5 g
X3277	C11-079	Start up kit for 1000 oxygen analyses 25/19 mm <ul style="list-style-type: none"> X1079V (1 each): Prepacked oxygen pyrolysis tube 25/19 mm (1pc) B2134 (1 each): Standards kit CHN/CHN-O, acetanilide, EDTA, urea, aspartic acid, 4 x 0.5 g D2009 (10 each): Silver capsules pressed 8 x 5 mm, pack of 100
X3278	C11-080	Consumable kit for 1000 analyses (CHNS w heteroatoms) 25/19 mm (CHNS) <ul style="list-style-type: none"> X1076V (1 each): Packed reactor 25/19 mm CHNS (1pc) D1008 (4 each): Tin capsules pressed standard weight 8 x 5 mm, pack of 250
X3279	C11-081	Consumable kit for 1000 oxygen analyses 25/19 mm <ul style="list-style-type: none"> X1079V (1 each): Prepacked oxygen pyrolysis tube 25/19 mm (1pc) D2009 (10 each): Silver capsules pressed 8 x 5 mm pack of 100

These products offer the same high quality and value for money that is typical of Elemental Microanalysis' entire product range. To place an order, call us, send an email to enquiries@microanalysis.co.uk, or register to view our product pricing and to purchase online. Elemental Microanalysis has over 45 years' experience manufacturing and supplying quality products with exceptional service and affordable pricing.

Elemental Microanalysis Ltd

1 Hameldown Road, Okehampton, EX20 1UB, UK

• t. +44(0)1837 54446 • www.elementalmicroanalysis.com

• enquiries@microanalysis.co.uk



Available from your local Dealer:

