

New Products for Elementar®

Our research and development, quality, product management and manufacturing teams collaboration has enabled Elemental Microanalysis to add these high quality, great value, compatible products to our product range.

A new range of compatible products for Elementar® analysers, including the inductor series have been introduced.

All Elemental products are available to order from our website.



Part description	Pack size	Part no.	OEM part no.	Compatible analysers
Combustion tube, 32 mm	1	C3073	100000245	Elementar® inductor CS cube/EL cube
Filter tube 42 mm	1	C3074	16.00-1081	Elementar® vario TOC cube / TOC select
Glass tube GL25, 140 mm	1	C3075	(100000465) glass part only	Elementar® inductor CS cube/EL cube/ONH cube
Glass tube GL25, 200 mm	1	C3076	(100000466) glass part only	Elementar® inductor EL cube/ONH cube
Glass tube GL25 140 mm	1	C3077	100000465	Elementar® inductor CS cube/EL cube/ONH cube
Glass tube GL25 200 mm	1	C3078	100000466	Elementar® inductor EL cube/ONH cube
Graphite crucible	500	C4648	100002030 ELCUP ONH	Elementar® inductor EL cube/ONH cube
O-ring, black 4 x1 mm	1	E1179	200011633 03002339	Elementar® inductor CS cube/EL cube/ONH cube
O-ring, green, 19.5 x 2.5 mm	1	E1180	100000341	Elementar® inductor CS cube/EL cube
O-ring, green, 24.5 x 3 mm	1	E1181	100000331	Elementar® inductor CS cube/EL cube
O-ring, white	1	E1182	100002419	Elementar® inductor EL cube/ONH cube
O-ring, 10.5 x 2mm	1	E1183	200009016	Elementar® inductor CS cube/EL cube/ONH cube

The use of Elementar® part numbers is intended for convenience only, and does not imply the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All Elementar® trademarks acknowledged.

Elemental Microanalysis Ltd

1 Hameldown Road, Okehampton, EX20 1UB, UK
 • t. +44(0)1837 54446 • www.elementalmicroanalysis.com
 • enquiries@microanalysis.co.uk



Available from your local Dealer:

