

## Now offering ICP-MS cones

Our **45 years** of expertise in offering high quality OEM alternative consumables to our customers is now extended to cover your ICP-MS needs! Elemental Microanalysis is now extended to cover your ICP-MS needs (cones and skimmers) that seamlessly support the needs of our customers who have trace elemental analysis capability laboratories.



**Micro-Skimmer Cone (V4001)**

As the leading independent supplier of elemental analysis consumables, we can offer a substantial savings in comparison to the prices charged by the original supplier of your instrument.

Ordering with Elemental Microanalysis gives you the assurance of quality while providing price savings for your lab.

OEM No.	EMAL No.	Description	Material	Suitable for the following Thermo Fisher Scientific® instruments			
				PQ2	VG ExCell	X-SERIES	iCapQ/Qnova
3201101	V4000	Micro-Skimmer Cone	Platinum	x	x	x	
3200860	V4001	Micro-Skimmer Cone	Nickel	x	x		
3004661	V4002	Sample Cone	Nickel	x	x	x	
3601289	V4003	Sample Cone	Platinum			x	x
3600812	V4004	Sample Cone	Nickel			x	x
3600811	V4005	Skimmer Cone	Nickel			x	
3601320	V4006	Skimmer Cone Adapter	Nickel			x	

OEM No.	EMAL No.	Description	Material	Suitable for the following Thermo Fisher Scientific® instruments		
				ELEMENT	ELEMENT 2	NEPTUNE
1044530	V4007	Sample Cone	Nickel	x	x	x
1067600	V4008	'H' Skimmer Cone	Nickel	x	x	x
1067500	V4009	Sample Cone	Platinum	x	x	x
1067501	V4010	Sample Cone (w/ Boron)	Platinum	x	x	x
1047460	V4011	'H' Skimmer Cone	Platinum	x	x	x
1047461	V4012	'H' Skimmer Cone (w/Boron)	Platinum	x	x	x
1142160	V4013	'X' Skimmer Cone	Nickel	x	x	x
1142150	V4014	'X' Skimmer Cone (Boron free)	Platinum	x	x	x
1260630	V4015	Jet Sample Cone	Nickel		x	x

# Guaranteed quality


Our confidence in our products allows us to offer you a total satisfaction guarantee: if our product does not perform at least as well as the brand of product you usually use, you can return it for a full refund.



**'H' Skimmer Cone (V4008)**



**Sample Cone (V4007)**



**Sample Cone (V4004)**



**Skimmer Cone (V4005)**

OEM No.	EMAL No.	Description	Material	Suitable for the following Thermo Fisher Scientific® instruments							
				VG ExCell	X-SERIES	iCapQ/ Qnova	ELEMENT	ELEMENT 2	NEPTUNE	Apix	
3200456	V4016	Plasmatrace Skimmer Cone	Nickel	x							
3201941	V4017	Skimmer Cone (7.2mm insert)	Platinum	x							
3201585	V4018	Skimmer Cone (0.4 Hole – Eclipse)	Nickel	x							
1067510   TS067016	V4019	Skimmer Cone (N)	Platinum				x	x	x		
1044540   TS067015	V4020	Skimmer Cone (N)	Nickel				x	x	x		
1044850   FS067017B	V4021	Skimmer Fitting Tool	Resin				x	x	x		
1142151	V4022	'X' Skimmer Cone (with Boron)	Platinum				x	x	x		
1184050	V4023	Sample Cone	Aluminium				x	x	x		
1184060	V4024	Skimmer Cone	Aluminium				x	x	x		
3600813	V4025	Skimmer Cone Adapter (w/4 screws)	Copper		x						
4600347	V4026	Xi Interface Kit			x						
3601219	V4029	Torch Shield	Silver		x						
3600814	V4030	Screws for Skimmer Cone Adaptor (pk/4)	Stainless Steel		x						
3601245	V4031	Skimmer Cone Adaptor	Copper Nickel		x						
1341430	V4032	Skimmer Cone (Hot/ Cold Plasma)	Platinum			x					
3202521	V4033	Sample Cone, Insert (6.35mm)	Platinum	x							
1600335	V4034	Sampler	Platinum		x						