

# Technical Information Bulletin

## Ion chromatography vials and standards

Elemental Microanalysis offers a range of high quality vials, caps and accessories for Thermo/Dionex® and Metrohm® Ion Chromatography auto-samplers at attractive prices.

The vials are manufactured from high grade polypropylene — the preferred material for low and consistent blank anion, cation, organic acid and Cr(VI) values. These products are available with plain, filtered and pierce caps and are fully compatible with the equipment mentioned. For details see the table on the following page. Available now for immediate delivery.



Elemental Microanalysis Ltd  
1 Hameldown Road Okehampton  
Devon EX20 1UB UK  
t. +44 (0) 1837 54446  
f. +44 (0) 1837 54544  
info@microanalysis.co.uk  
[www.elementalmicroanalysis.com](http://www.elementalmicroanalysis.com)



4260



® - All trademarks acknowledged. Please note that manufacturer part numbers are given for reference only and do not imply that the products are of Thermo/Dionex or Metrohm origin. We offer high quality alternatives of equivalent performance.

Available from your local Dealer:

# Vials for Dionex & Metrohm IC

OEM Pt No Dionex/Metrohm	Pt No	Volume	Description	Application	Pack
<b><u>DIONEX VIALS</u></b>					
<b><u>5ml Vials</u></b>					
038141	F7000	5ml	5ml Vial with Cap (with filter) for Dionex® AS40/ASDV	Dionex AS40/ASDV	250
039532	F7001	5ml	5ml Vial with Cap (without filter) for Dionex® AS40/ASDV	Dionex AS40/ASDV	250
038008	F7005	5ml	5ml Vial (Vial Only) for Dionex® AS40/ASDV	Dionex AS40/ASDV	250
038009	F7006		Cap (with filter) for Dionex® AS40/ASDV 5ml Vials	Dionex AS40/ASDV	250
039528	F7007		Cap (without filter) for Dionex® AS40/ASDV 5ml Vials	Dionex AS40/ASDV	250
<b><u>0.5ml Vials</u></b>					
038142	F7002	0.5ml	0.5ml Vial with Cap (with filter) for Dionex® AS40/ASDV	Dionex AS40/ASDV	250
	F7003	0.5ml	0.5ml Vial with Cap (without filter) for Dionex® AS40/ASDV	Dionex AS40/ASDV	250
038010	F7009	0.5ml	0.5ml Vial (Vial Only) for Dionex® AS40/ASDV	Dionex AS40/ASDV	250
038011	F7010		Cap (with filter) for Dionex® AS40/ASDV 0.5ml Vials	Dionex AS40/ASDV	250
079797	F7011		Cap (without filter) for Dionex® AS40/ASDV 0.5ml Vials	Dionex AS40/ASDV	250
<b><u>1.5ml Vials</u></b>					
079812	F7004	1.5ml	Vial Kit, 1.5ml Polypropylene with Caps and Septa, for Dionex® ASAP/AS/AS50	Dionex AS50/AS/ASAP	100
<b><u>7–10ml Vials</u></b>					
Alternative for 055058	F7012	7ml	7ml Vial with Pierce Cap for Dionex® AS/ASAP/AS50	Dionex AS50/AS/ASAP	250
<b><u>METROHM VIALS</u></b>					
6.2743.050 (2000 pcs)	F7015 (1000)	11ml	Sample tubes 11ml for IC Sample Processors and VA Autosampler from Metrohm	Metrohm IC	1000
6.2743.057	F7016	11ml	Sample tubes 11ml for IC Sample Processors and VA Autosampler from Metrohm	Metrohm IC	200
6.2743.070 (2000 pcs)	F7017 (1000)		Stopper with perforation for 11ml sample tube	Metrohm IC	1000
6.2743.077	F7018		Stopper with perforation for 11ml sample tube	Metrohm IC	200