


# *New product release*

## **Replacement Furnace for Thermo Flash EA1112**



A quality alternative to the part offered by the original manufacturer. Adaptors are supplied to configure the furnace for 18mm or 25mm reactors.

This product is equivalent to Thermo part numbers 354 60018 and 354 60021.

### **Ordering details:**

<i>Cat No</i>	<i>Description</i>	<i>Pack</i>
E4410	Replacement Furnace Flash EA	1

Elemental Microanalysis  
1 Hameldown Road  
Okehampton  
Devon, EX20 1UB, UK  
+44 (0)1837 54446  
[info@microanalysis.co.uk](mailto:info@microanalysis.co.uk)  
[www.microanalysis.co.uk](http://www.microanalysis.co.uk)

Or your local dealer:

Apr 08

All trade marks acknowledged. The use of manufacturer's part numbers is for convenience only and does not imply that the products are of the manufacturer's origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives.

**Replacement Furnace for Flash EA1112**