


New product release

Large Tin Foil Cones



D9002 Large tin foil cones

D9001 standard tin foil cones

A larger size cone formed from a 50mm suitable for large volume samples. This product is equivalent to Leco® part numbers 502-397 and 502-397-400.

Ordering details:

Cat No	Description	Pack
D9002	Large Tin foil Cones	100
D9003	Large Tin foil Cones	400

Elemental Microanalysis
 1 Hameldown Road
 Okehampton
 Devon, EX20 1UB, UK
 +44 (0)1837 54446
info@microanalysis.co.uk
www.microanalysis.co.uk

Or your local dealer:

Apr 08

All trade marks acknowledged. The use of manufacturer's part numbers is for convenience only and does not imply that the products are of the manufacturer's origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives.

Large Tin foil Cones