

LITERATURE REQUEST

Elemental Microanalysis can supply product sheets tailored specifically to your analyser, product catalogues, and various other literature to help you find the products you need. This literature is available in printed or digital format – just tick the boxes, fill out your details and return the form to us or your local dealer. The latest editions of our Catalogues are also available to download from our website: www.elementalmicroanalysis.com

Product sheets available, by analyser:

Inorganic *(Please Tick)*

Bruker G4 Icarus CS	
Bruker G8 Galileo O, N, H	
Eltra CS2000	
Eltra CS800	
Eltra OH900	
Eltra ON900	
Horiba EMGA-620/621	
Horiba EMGA-800/900	
Horiba EMGA-821/921	
Horiba EMIA-320V2/920V2	
Horiba EMIA820V/320V/220V	
Leco 836/736 ONH/OH/NH/H	
Leco CS125/255	
Leco CS200/C200/S200/WC200	
Leco CS230	
Leco CS244/344/IR12/232	
Leco CS300	
Leco CS400/C400/S400	
Leco CS444/IR412/432	
Leco CS600	
Leco CS744/844	
Leco DH102/103/RH101	
Leco DH603	
Leco RH402	
Leco RH404	
Leco RHEN602	
Leco RO478	
Leco TC136/336/TN114/RO116	
Leco TC300/RO300/TN300	
Leco TC436/TN414/RO416	
Leco TC500/600	
Leco TCH600	
Leco WR112	
Strohlein/Juwe NH-mat 243	

Micro/CHN(OS) *(Please Tick)*

Costech ECS4010	
Elementar CHN Rapid	
Elementar Pyro Cube	
Elementar Vario Micro Cube	
Elementar Vario EL	
Elementar Vario EL Cube	
Elementar Vario Isotope Cube	
Europa Roboprep/ANCA/Integra	
Eurovector EA3000	
Eurovector EA3100	
Eurovector EA-3028HT	
Exeter CE240/440	
Hekatech HT-O	
Leco 900/932	
Leco Truspec Micro	
Perkin Elmer 2400	
Perkin Elmer 240ABC	
Picarro CM Model 02	
SerCon GSL/Integra	
Thermo EA1106CHNSO	
Thermo EA1108CHNSO	
Thermo EA1110CHNSO	
Thermo EA1112CHNSO/FLASH2000	
Thermo Flash 2000 IRMS	
Thermo Flash HT/HT Plus	
Thermo Flash Smart CHNOS	
Thermo NA1500NCS2	
Thermo TC/EA	
Thermo NA2500NCS	

Macro CHN/CNS *(Please Tick)*

Elementar rapid CS cube	
Elementar Vario Macro	
Elementar Vario Macro cube	
Elementar Vario Max CHN/CNS	
Elementar Vario Max Cube N/CN	
Elementar Vario Max N/CN	
Leco CHN1000/CNS1000	
Leco CHN2000	
Leco CHN600/800	
Leco CNS/NS/CN2000	
Leco Truspec CHN/CN	
Skalar Primacs	

N/Protein *(Please Tick)*

Buchi DuMaster D-480	
CKIC TCN/CHN 2200	
Costech Turbo N	
Elementar Rapid N®	
Elementar Rapid N Cube®	
Elementar Rapid N exceed®	
Elementar Rapid MAX N exceed®	
Gerhardt Dumatherm®	
Leco 828 series	
Leco FP2000	
Leco FP228	
Leco FP328	
Leco FP428	
Leco FP528	
Leco FP628/CHN628	
Leco FP/CN928 series	
Leco TruMac	
Leco Truspec N	
Perkin Elmer 2410	
Thermo Flash 4000	
Thermo FlashSmart CHNOS/N-Protein	
Thermo NA1112Protein/FLASH2000	
Thermo NA2000/2100	
Velp NDA 701 Dumas Nitrogen Analyser	

Organic Carbon/Sulphur

(Please Tick)

Eltra CHS-580/CHS-580A	
Eltra CS-500	
Eltra Helios	
Eltra SurfaceC-800	
Leco SC632	
Leco SC832 series	
Leco TGA601/701	
Leco RC612	
Leco SC144	
Leco SC32/132/232/CR12	
Leco SC432/444/CR412	
Leco RC412	

All trademarks acknowledged. Use of the OEM part number is for reference only and does not imply that the products are of OEM origin. All Elemental Products are guaranteed to be of high quality and suitable for use in the stated application.

Product sheets (continued)

Trace (Please Tick)

Analytik Jena Multi NC2100-3100	
Analytik Jena Multi EA2000	
Analytik Jena Multi EA3000	
Analytik Jena Multi EA4000	
Analytik Jena Multi EA5000	
Analytik Jena Multi X2000	
Analytik Jena Multi X2500	
Analytik Jena Pharma TOC	
Antek 7000 from PAC	
Antek 9000 from PAC	
Antek MultiTek analyser from PAC	
Dohrmann DC	
Dohrmann DN	
Dohrmann MCTS	
Dohrmann DX	
Elementar LiquiTOC	
Elementar trace SN cube	
Elementar Vario TOC Cube	
Elementar Vario Trace S	
Hach IL550 TOC-TN analyser	
Metrohm 930 Combustion IC PP	
Shimadzu TOC	
TE Instruments	
Thermo Euroglas All Instruments	
TSHR Instruments	

Other literature available:

Catalogues (Please Tick)

Capsules	
Inorganic	
IRMS Supplement	
Laboratory Analysis Service	
Macro CHN/CNS (including configuration guide)	
Micro/CHN (OS)	
Nitrogen / Protein	
Organic Carbon / Sulphur	

Reference materials/Standards for EA

Pure chemical standards	
Natural materials	
Inorganic materials	
Isotope reference materials	

Ion Chromatography (Please Tick)

Dionex AS40/ADSV/ASAP/AS50/AS IC vials	
Metrohm Peak IC Autosampler IC vials	

Can't find what you are looking for? Let us see if we can help: info@microanalysis.co.uk or +44 (0)183754446

What format would you like your literature in?
(Please circle)

Print

Digital

Name

Organisation

Address

Tel

Email

Post, Fax or Email to:
Elemental Microanalysis Ltd
1 Hameldown Road
Okehampton, Devon
EX20 1UB UK
F: +44 (0)1837 54544
info@microanalysis.co.uk

Or Your Local Dealer: