

New product release

Economy Aluminium Moisture Determination Pans



D5017

An economical alternative to our original pans, D5014, these new Aluminium Moisture Determination Pans are manufactured from high quality material and still retain the same measurements 100mm diameter and 7mm high but cost a good deal less. Suitable for most types of moisture determination techniques.

For even bigger savings try our multipack – D5018 see below for details.

Ordering details:

Cat No	Description	Pack
D5017	Economy Aluminium Moisture Pans 100 dia x 7mm high	80
D5018	Economy Aluminium Moisture Pans 100 dia x 7mm high	10x80

Elemental Microanalysis
 1 Hameldown Road
 Okehampton
 Devon, EX20 1UB, UK
 +44 (0)1837 54446
info@microanalysis.co.uk
www.microanalysis.co.uk

Or your local dealer:



Jul 2006

All trade marks acknowledged. The use of manufacturer's part numbers is for convenience only and does not imply that the products are of the manufacturer's origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives.

Aluminium Moisture Determination Pans