

Encapsulation for Elementar®

Elemental Microanalysis supplies customers all over the world with encapsulation products compatible with Elementar instruments.

Available to order from our ISO 9001 factory, here are some of our customer favourites.



Tin foil cones standard weight

Elementar	EMAL
05 000 429	D9001
05 002 797	D9010



Smooth wall tin capsules round base 20 x 8mm

Elementar	EMAL
05 000 311	D4009



Tin foil squares standard weight 37 x 37mm and 50 x 50 mm

Elementar	EMAL
12.01-0032	D1116
12.01-0033	D1115 D1118
34.00-0052	D1119



Tin capsules pressed standard weight 5 x 3.5mm

Elementar	EMAL
05 003 395	D1002



Tin boats 15 x 8 x 8mm

Elementar	EMAL
22 137 420	D5032



Silver capsules pressed 8 x 5mm

Elementar	EMAL
05 003 397 05 001 540	D2008




Tin boats 12 x 4 x 4mm

Elementar	EMAL
22 137 418	D5130



Smooth wall nickel capsule flat base 14.5mm x 7mm

Elementar	EMAL
ELCAPS OHN	D3509



Silver boats 12 x 6 x 6mm

Elementar	EMAL
22 133 213	D5034

The use of Elementar® part numbers is intended for convenience only, and does not imply the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All Elementar® trademarks acknowledged.

Encapsulation for Elementar®

Elementar part no.	Material	Description	Pack size	EMAL part no.
03 951 620	Tin	Smooth wall capsules flat base 6 x 3mm	250	D4131
05 000 311	Silver	Smooth wall capsules round base 20 x 8mm	100	D4009
05 000 312	Tin	Smooth wall capsules flat base 14 x 6.75mm	100	D4120
05 000 429 05 002 797	Tin	Foil cones standard weight	100 1500	D9001 D9010
05 000 481	Tin	Smooth wall capsules flat base 11 x 4.75mm	250	D4121
05 000 816	Silver	Capsules smooth wall 5.5 x 3.5mm	250	D4102
05 001 727	Tin	Smooth wall capsules flat base 5.5 x 3.5mm	100	D4033
05 003 394	Tin	Capsules pressed standard weight 8 x 5mm	250	D1008
05 003 395	Tin	Capsules pressed standard weight 5 x 3.5mm	500	D1002
05 003 396	Silver	Capsules pressed 5 x 3.5mm	500	D2002
05 003 397 05 001 540	Silver	Capsules pressed 8 x 5mm	250	D2008
05 003 398	Silver	Capsules pressed 11 x 8mm	100	D2103
05 003 400	Tin	Capsules pressed standard weight 11 x 8mm	100	D1103
12.01-0032	Tin	Foil squares standard weight 37 x 37mm	400	D1116
12.01-0033 12.01-0033 34.00-0052	Tin	Foil squares standard weight 50 x 50mm	100 500 1500	D1115 D1118 D1119
22 133 212	Silver	Weighing boats 12 x 6 x 6mm	100	D5026
22 133 213	Silver	Weighing boats 15 x 8 x 8mm	70	D5034
22 137 418	Tin	Boats 12 x 4 x 4mm	1000	D5130
22 137 419	Tin	Boats 12 x 6 x 6mm	500	D5131
22 137 420	Tin	Boats 15 x 8 x 8mm	250	D5032
25.00-0074	Silver	Squares 35 x 35mm	100	D2084
ELCAPS OHN	Nickel	Smooth wall capsule flat bottom 14.5mm x 7mm	100	D3509

The use of Elementar® part numbers is intended for convenience only, and does not imply the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All Elementar® trademarks acknowledged.

Elemental Microanalysis Ltd

1 Hameldown Road, Okehampton, EX20 1UB, UK
t. +44(0)1837 54446 • www.elementalmicroanalysis.com •
enquiries@microanalysis.co.uk



Available from your local Dealer:

