



Product Sheet ~ Leco Truspec N

~ Other sample encapsulation

Cat #	Description	OEM #
D9005	Cel Capsules Medium 502-382 pack of 400	502-382
D9004	Cel Capsules Small 502-338 pack of 400	502-338

~ Tin capsules for liquids

Cat #	Description	OEM #
D8012	Evaporating Foil Cups Tin Capacity 5ml 28 x 24mm diameter pack of 48	UK-338
D4003	Smooth Wall Tin Capsules Round Base 18 x 5mm pack of 100	501-059
D4015	Smooth Wall Tin Capsules Round Base 18 x 6mm 502-040-100 pack of 1000	502-040-100
D4007	Smooth Wall Tin Capsules Round Base 18 x 6mm pack of 100	502-040
D4009	Smooth Wall Tin Capsules Round Base 20 x 8mm pack of 100	502-167

~ Tin cones squares and discs

Cat #	Description	OEM #
D9002	Large Tin Foil Cones Standard Weight 502-397 pack of 100	502-397
D9003	Large Tin Foil Cones Standard Weight 502-397-400 pack of 400	502-397-400
D9001	Tin Foil Cones Standard Weight 502-186 pack of 100	502-186 502-186-200
D9010	Tin Foil Cones Standard Weight 502-186B pack of 1500	502-186-100

~ Tools, trays etc

Cat #	Description	OEM #
E2016	Forceps Stainless Steel Curved Points	604-378
E2017	Forceps Stainless Steel Straight Points	502-028
E2462	Spatula 501-614	501-614

Chemical Reagents

Cat #	Description	OEM #
B1163	Aluminium Oxide Pellets 502-188 200gm	502-188
B1107	Copper Turnings Fine Long Strands Reduced and Annealed 501-621/502-656 200gm	501-621 502-258 502-295 502-656
B1015M	Copper Wires Fine Wires Reduced 4 x 0.5mm 100gm	502-189 502-304 502-705 502-844
B1212M	Copper Wires Fine Wires Reduced 4 x 0.5mm 1kg	-
B1083M	Copper Wires Pound Pack Fine Wires Reduced 4 x 0.5mm 454gm	502-304-500
B1161	EMADrone Magnesium Perchlorate 501-171 454gm	501-171 501-171-HAZ
B1162	EMASorb A Granular 0.8 to 1.6mm 14 to 25 mesh 502-174 500gm	502-174 502-174-HAZ

B1164	EMASorb B Granular 1.6 to 3mm 7 to 14 mesh 502-176 500gm	502-176
B1191	Furnace Reagent 501-609 100gm	501-609
B1166	Glass Wool Borosilicate Fibres 501-081 454gm	501-081
B1167	Glass Wool Borosilicate Fibres 763-265 227gm	763-265
B1266	Glass Wool Rope Borosilicate Fibres 501-081 454g	501-081
B1267	Glass Wool Rope Borosilicate Fibres 763-265 227g	763-265
B1195	Iron Chip Accelerator 501-077 908gm	501-077
B1165	Iron Wool Grade #2 502-310 450gm	502-310
B1190	Magnesium Oxide Beads 502-359 200gm	502-359
B1061	Nitrogen Catalyst Pellets 502-049 50gm	502-049
B1192	Silica (Quartz) Wool Strips	608-379

Crucibles and Boats

Cat #	Description	OEM #
C4038	Ceramic Boats Unglazed 529-203 pack of 50	529-203
C4071	LARGE CRUCIBLE FOR LIQUID 617-605 Pack of 10	617-605

Natural Reference Materials

Cat #	Description	OEM #
B2164	Algae (Bladderwrack) Standard 30gm See Cert.117020 Expiry 21-July-21	-
B2162	Algae (Spirulina) Standard 30g See Certificate 290349 Expiry 28-Sep-20	-
B2166	Birch leaf Standard 30g See Certificate 304941 Expiry 22-Mar-22	-
B2160	Chitin Standard 30g See Certificate 293513 Expiry 15-Nov- 20	-
B2320	Coal CHNS Standard (S~0.41%) RM See Certificate 751219 for Values 50gm	502-383
B2321	Coal CHNS Standard (S~1.12%) RM See Certificate 776419 for Values 50gm	502-384
B2322	Coal CHNS Standard (S~2.27%) RM See Certificate 781918 for Values 50g	502-385
B2324	Coal CHNS Standard S=5.39 RM See Certificate 782118 for Values 50gm	502-388
B2168	Coconut shell Standard 30g See Certificate 208484 Expiry 23- Sep- 23	-
B2154	Protein Standard OAS See Certificate 199622 Exp.11-May-2023(Certified for Nitrogen only) 30g	-
B2158	Sorghum Flour Standard OAS 30gm See Certificate 114856 Expiry 17-March 2023	-
B2156	Wheat Flour Standard OAS 30g See Certificate 280967 Exp 11-May-23	-

O-rings and operational accessories

Cat #	Description	OEM #
E2203	Cellulosic Tubes (Sampling Tubes) 502-210 pack of 500	502-210
E2794	C-Flex Tubing Repl Kit 625-602-172	625-602-172 625-602-470
E2154	Cleaning Brush 601-981	601-981
E2165	Cooling Plate Stainless Steel 601-387	601-387
E2199	Lance Assembly Complete 616-147	616-147
E2197	Lance Holder Block 783-275	783-275

Product Sheet ~ Leco Truspec N

E2198	Lance Holder Cover 783-276	783-276
E2379	Mesh screen filter 617-441	617-441
E1125	Neoprene Gasket Ballast Seal 778-120, 160mm x 150mm x 3mm	778-120
E1344	O Ring 601-684, 31.4mm x 2.6mm	601-684
E1102	O Ring 601-939 4.7mm x 1.4mm	601-939
E1345	O Ring 617-136, 5.2mm x 2.6mm	617-136
E1311	O Ring 771-407, 17.1mm x 2.6mm	771-407
E1342	O Ring 774-425, 41.0mm x 2.6mm	774-425
E1358	O Ring 619-139, 25.1mm x 2.6mm	619-139
E1359	O Ring 619-192, 56.9mm x 1.8mm	619-192
E1360	O ring 702-557, 29.8mm x 2.6mm	702-557
E1110	O Ring 772-520, 6.0mm x 2.6mm	772-520
E1116	O Ring Ballast Tank 778-116, 145.7mm x 3.5mm	778-116
E1111	O Ring Catalyst Heater Assembly 772-738, 6.1mm x 1.8mm	772-738
E1101	O Ring Catalyst Heater Tube 601-442 14.0mm x 2.6mm	601-442
E1107	O Ring Condensing Block 763-379 44.2mm x 1.8mm	763-379
E1109	O Ring Doser Valve Assembly 772-269 pack of 10, 2.9mm x 1.8mm	772-269
E1133	O Ring Filter Tube 765-976, 15.5mm x 3.5mm	765-976
E1114	O Ring Sample Drop Block 773-913, 4.5mm x 1.8mm	773-913
E1140	O-Ring Combustion Tube 616-138, 36.1mm x 3.5mm	616-138
E2464	O-Ring Installation Tube 603-350	603-350
E2369	Particle Filter 616-504	616-504
E2150	Pinch Valve Tubing 601-69 ft	601-691
E2163	Screw Nylon 6-32 x 1 Slotted Countersunk Head 190-865	190-865
E2201	Screws For Lance Assembly 190-423 pack of 4	190-423
E2235	Secondary (screen) Filter 775-306	775-306
E2158	Silicone Grease 501-241 100gm	501-241
E2039	Snoop Leak Detection Fluid 236ml/8 fluid oz	502-213
E2167	Spring Dosing Valve Assembly 602-672	602-672
E2182	Tube Fitting Brass 10/32 to 1/16 ID tube 773-661	773-661
E2181	Tube Fitting White Nylon 1/8 NPT to 3mm ID tube 602-605	602-605

Prepacked reactors and consumables kits

Cat #	Description	OEM #
X5151	Consumable Kit 1000 Analysis for TruEssential (Small Tin Foils) 502-445-HAZ	502-445-HAZ
X5150	Consumable Kit 5000 Analysis for Truspec (Small Tin Foils) 502-621- HAZ	502-621-HAZ
X3000	Prepacked Catalyst Heater/Reduction Tube Leco Transparent Silica	-
X3015	Prepacked Reaction Tube Leco Truspec	-

Pure Chemical Standards

Cat #	Description	OEM #
B2144	Ethylenediaminetetra-acetic Acid (EDTA) 502-092B 500gm See Certificate 335390 Expiry 14-Jul-24	502-092-250
B2137	Ethylenediaminetetra-acetic Acid (EDTA) OAS 100gm See Cert 335390 Expiry 14-Jul-24	-
B2115	Ethylenediaminetetra-acetic Acid (EDTA) OAS 502-092 50gm See Certificate 335390 Expiry 14-Jul-24	502-092
		502-812
B2121	Sucrose OAS See Certificate 201037 Expiry 04-May-22 501-046 25gm	501-046
		501-441

Quartz, Glassware and Furnace tubes

Cat #	Description	OEM #
C4032	Bubble/Filter/Reagent Tube Borosilicate Glass 771-030	771-030
C4023	Catalyst Heater Tube Transparent Silica 601-437	601-437
C4123	Catalyst Heater Tube Transparent Silica 601-437B pack of 5	-
C4338	Combustion Tube 619-065	619-065
C4027	Filling Funnel Borosilicate Glass 502-023	502-023
C4342	Filter tube 618-809	618-809
C4037	Furnace Filter Tube Borosilicate Glass 601-390	601-390
C4060	Glass Scoop 503-032	503-032
C4058	Honeycomb Ceramic 780-899	780-899
C4051	Oxygen Lance Long (Lance length 5.25") 616-146	616-146
C4066	Oxygen Lance Short (lance length 4.75) 614-962	614-962
C4065	Porous Crucible (Pack of 10) 614-961-110	614-961 614-961-110
C4067	Porous Crucible (Pack of 5) (1/2 pack of 614-961-110)	-
C4332	Reagent Tube 617-548/633-103-225	617-548 633-103-225
C4341	Reagent tube 619-268	619-268
C4004	Reagent Tube Borosilicate Glass 775-601	775-601
C4339	Reduction Tube 619-154	619-154

Quality: All products are manufactured in our ISO9001:2015 approved facility and are covered by our 'Total Satisfaction Guarantee'.

Can't see what you need? We offer thousands of products for all elemental analysers. If you can't see what you need here the chances are we will have the product available. Just contact our friendly sales team who will track down what you need.

OEM Part numbers: The use of OEM part numbers is intended for convenience only and does not imply that the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and suitable for use in the stated instrument. All trade marks acknowledged.