

# New catalysts, crucibles and furnace tubes



Our research and development, quality, product management and manufacturing teams collaboration has enabled Elemental Microanalysis to add these high quality, great value, compatible products to our product range.

Featuring in our new release of the Micro catalogue the **C1038** joins our growing range of products compatible with the Sercon isoEarth<sup>+</sup> analyser.

Three new products for the Elementar<sup>®</sup> Isoprime TraceGas module have been added; **C1114** is the GC/IRMS ceramic furnace tube, **X1508** is a prepacked version of this furnace tube; it is supplied with the copper, platinum and nickel chromium wire expertly inserted into the narrow tube.

**Prefer the challenge of packing the GC/IRMS ceramic furnace tube? D9103** (Nickel/chromium 80/20 wire 0.125 mm) is now available to purchase alongside existing products **D9100** (0.125 mm platinum wire) and **D9101** (0.125 mm copper wire).

Customer request has instigated the release of **D4603-1000**, a 1000 unit pack size of cylindrical graphite crucibles compatible with LECO<sup>®</sup> RH404, LECO<sup>®</sup> RHEN602.

Combustion tube **C4702** (3 200 607 909 equivalent) is Elemental's first compatible product for the Horiba<sup>®</sup> EMIA Expert / EMIA Pro Inorganic analyser series. All these products and more are available to order from our website.



D9103



C4603-1000



C4702

C1038

C1114 and X1508

Part description	Pack size	Part no.	OEM part no.	Compatible analysers
Quartz reduction tube, 6 mm OD, 4 mm ID, 370 mm	1	C1038	SC9131	Sercon <sup>®</sup> isoEarth <sup>+</sup>
GC/IRMS furnace reaction tube ceramic, 450 mm x 3 mm OD x 0.5 mm ID	1	C1114	SM679095AD2	Elementar <sup>®</sup> Isoprime TraceGas
Graphite crucible	1000	C4603-1000	764-330	LECO <sup>®</sup> RH404, LECO <sup>®</sup> RHEN602
Combustion tube for Horiba <sup>®</sup>	1	C4702	3 200 607 909	Horiba <sup>®</sup> EMIA Expert / EMIA Pro
Nickel/chromium 80/20 wire 0.125 mm	10 m	D9103	ST104051	Elementar <sup>®</sup> Isoprime TraceGas
Prepacked C1114, 260 mm copper / platinum / Nichrome wire (260 mm)	1	X1508	None Available	Elementar <sup>®</sup> Isoprime TraceGas

The use of OEM part numbers is intended for convenience only, and does not imply the product of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All OEM trademarks acknowledged.

## Elemental Microanalysis Ltd

1 Hameldown Road, Okehampton, EX20 1UB, UK  
 • t. +44(0)1837 54446 • [www.elementalmicroanalysis.com](http://www.elementalmicroanalysis.com)  
 • [enquiries@microanalysis.co.uk](mailto:enquiries@microanalysis.co.uk)



Available from your local Dealer:

