

Dealer support



Exhibition items available to dealers

Are you looking for a simple and professional way to exhibit our product range? Elemental Microanalysis can help you stand out at your trade show, exhibition, conference or when meeting customers.

Items can be loaned or purchased for display use, offering you an opportunity to physically demonstrate Elemental's product quality and range in a smart, yet compact format.


Marketing and promotional documents can be translated on request (*Elemental Microanalysis does not supply translations*)

To request items for your exhibition, please contact marketing@microanalysis.co.uk



Item	Part number	Purchase available	Loan available
Glass display stand 1 with lighting. With tubes C3040, C8060, C4338, C4054, C4701, C4055	Z6002	✓	✓
Glass display stand 2 with lighting. With tubes C1039, C1069, C1074, C2015, C2016, C3009	Z6003	✓	✓
Reagent vials display (labelled)	Z6004	✓	✓
Display box 1: encapsulation (discs/squares/capsules)	Z2039	✓	✓
Display box 2: encapsulation (boats/capsules/cones)	Z2040	✓	✓
Display Box 3: encapsulation (ceramic/graphite/crucibles)	Z2041	✓	✓
Quartz wool discs display box	Z6006	✓	✓
Dionex® compatible vials display box	Z6007	✓	✓
ICP-MS cones Agilent Technologies® compatible	Z2047	✓	✓
ICP-MS cones NU Instruments compatible	Z2048	✓	✓
ICP-MS cones PerkinElmer® compatible	Z2049	✓	✓
ICP-MS cones Thermo Fisher Scientific® compatible	Z2050	✓	✓





Thanks for joining our dealer family!

Support materials




LinkedIn posts




E-mail updates




Product sheets



Dealer Sharepoint



Blog posts



Account page